


Baxter

AK 98

DIALYSIS MACHINE



1. Saline bag
2. Roller clamp or clamp on infusion line
3. Infusion line
4. Connector
5. Arterial blood line
6. Arterial blood line clamp
7. Arterial pressure transducer line
8. Arterial pressure transducer connector
9. Arterial infusion line
10. Heparin line
11. Syringe pump piston
12. Blood pump
13. Dialyzer
14. Venous blood line
15. Venous infusion lines
16. Venous drip chamber
17. Venous blood line clamp
18. Priming detector
19. Venous pressure transducer line
20. Venous pressure transducer connector
21. Dialysis fluid tubes
22. Dialysis fluid tube connectors
23. Waste bag

AK 98

Blood Circuit

Machine Introduction SW 3.XX

Important notice: The information in this document does not relieve any user from his or her duty to carefully read the full text of the Operator's Manual for the **AK 98** dialysis machine.

Baxter & AK 98 are trademarks of Baxter International Inc., or its subsidiaries.